
ingo-man® plus

Seifen- und Desinfektionsmittelspender
Dispenser for soap and disinfectant
Distributeur de savon et de désinfectant

Wichtig! Important!

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The compatibility of the soap/disinfectant material used in this dispenser must be checked prior use, in order to guarantee fault free functionality. OPHARDT Hygiene assumes no liability for any damages or costs which arise due to the incompatibility of filling material and dispenser systems.

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E Important!

Please carefully read the operating instructions and save them for reference. The proper function of this dispenser system can only be guaranteed if the procedures in these instructions are followed.

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
E Technical Drawing
Technical modifications reserved.

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


E 1. Preparation

When using as a soap dispenser, install the dispenser as close as possible to the sink to prevent wet hands from dripping onto the floor.

The dispenser is mounted to a vertical surface using 4 screws. An optional 2 screws can be installed to permanently secure the dispenser to a vertical surface, if required.

Required space above the dispenser: 50-70 mm for the lever, and for dispenser removal and placement, as well as an additional 100 mm for operation. The required room below the dispenser, necessary for bottle removal, is about 300 mm.

Optimal mounting height for use: 130-135 cm, measured from the floor to the upper edge of the dispenser housing.

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E 2. Mounting of the dispenser

First remove the bottle: to do this tilt the bottom of the bottle forwards and then pull the bottle out of the dispenser.

Remove the dispenser from the wall plate: press the latch on the back wall of the dispenser and slide the dispenser upwards, off the plate. With the enclosed, or other suitable attachment materials, attach the wall plate to the wall. Slide the dispenser onto the wall plate from the top until the latch clicks into place (see image 2).

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


⚠ Wall plates already in use from the ingo-man® classic/plus product family may continue to be used.
If mounting is necessary, please test wall for suitability before mounting.

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E 3. Flexible Sealing Cap

The included flexible sealing cap ensures optimal sealing for various bottle sizes and prevents the loss of alcohol from the disinfectant solution. When fitting a pump to the dispenser, always ensure that the bottle sealing cap is pulled down the pump cylinder at least 1 cm.

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E 4. Bottle Placement / Exchange

Insert the filled bottle into the dispenser from the bottom, until snaps into place.

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E 5. Initial Operation

For first use hold your hand under the nozzle and operate the lever several times until solution is dispensed.

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


E 6. Dosage

The dispenser is pre-set to dispense the maximum dosage of either 3.0 ml or1.5 ml/stroke, depending on the pump installed. To change the dosage, first remove the pump from the dispenser (see step 7).

⚠ **Important!** Do not take apart the housing! Use only original replacement parts issued from the manufacturer. Only for the dosage of suitable solutions such as soaps, lotions and disinfectants. Caution notice! The use of unsuitable solutions can influence the function of the dispenser and may cause the solution to spray into eyes.

The dosage can be changed by turning the adjustable wheel. On first use, the wheel may be difficult to turn; if necessary, use a pair of flat-nose pliers. The set dosage can be read from the wheel, when looking at it from the front, right side. The upper indication of the dosage volume is relevant (see image 6).

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

E 7. Pump Removal / Exchange


Gently pull face shield forward and then lift upwards. Next, press down on the pump lock with a finger, and pull the pump forwards (see images 7.1-7.3).

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
E 8. Cleaning / Sterilization

Remove the dispenser from the wall plate: to do this, press the latch on the back wall of the dispenser and slide the dispenser upwards, off the plate.

The ingo-man® plus stainless steel pumps are autoclavable with superheated steam (134 °C, 2 bar, max. 20 min dwell time).

⚠ **Important!** The dispenser, wall plate, bottle and plastic pumps (K pump) are not autoclavable! The dispenser is not dishwasher proof. For pump removal, see step 7. After autoclaving the stainless steel pump, allow it to cool for at least an hour prior to next use.

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E 9. Removal of Clip

If present, first remove the locking plate.

Remove the bottle: tilt the bottom of the bottle forwards and lift the bottle out. Remove the dispenser from the wall plate: to do this, press the latch on the back wall of the dispenser and slide the dispenser upwards, off the plate.

9.1 Turn the dispenser on its head, so you look into the housing through the bottom (see image 9.1).

9.2 Using your thumb, press the spring-loaded metal pin, located on the upper, back-side of the dispenser, in the direction of the pump. At the same time, use your other hand to carefully unhinge and remove the face shield from the front side of the dispenser.

9.3 Carefully press on the clip from the inside outwards, to release it from its socket (see image 9.2).

Assemble in reverse order (see image 9.3).

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E 10. Information for General Care


The Robert Koch-Institute recommends: "For both dispenser types (note: hand disinfectant and soap dispenser), the fixed outer and inner parts must be able to be reprocessed. This is also possible by machine. […] Visible contamination of the outer housing and all accessible parts without further manipulation as well as soilings at the outlet must be removed by wipe disinfection. Daily wipe disinfection is recommended for operating levers. The scope and frequency of checks and the reprocessing of dispensers should be specified in the internal hygiene plan, e.g. every six months." [Source: Robert Koch Institute, recommendations for hand hygiene, 2016]

In order to guarantee proper functioning of the dispenser regular manual cleaning or processing in cleaning-disinfection devices is necessary.
Manual cleaning: Regularly remove adherent product residual. Cleaning of the dispenser housing and the pump with hot water. Rinsing the pump with hot water. Disinfection of shroud is possible with appropriate surface disinfectants.
Processing in cleaning-disinfection devices: choose program „Dispenser cleaning" and use pH-neutral cleaning solutions. Important: Cleaning of plastic pumps up to max. 65 °C!
Processing in autoklave: the dispenser and stainless steel pump are autoclavable with superheated steam (134 °C, 2 bar, max. 20 min dwell time).
Plastic wall plate, bottle, plastic pumps (K pumps) are not autoclavable.

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E 11. Locking plate (optional)

11.1 Insert the top of the locking plate into the slot under the face shield. The key should be in open position. Push the bottom of the locking plate up and forward to close it in place. Turn the key clockwise until the key is vertical to lock the dispenser. Retract the key from the lock

11.2 To open the locking plate insert the key and turn fully anticlockwise, by applying upward pressure on the locking plate, pull forward.

11.3 After the locking plate is disengaged from a backplate pull it down to detach it from the dispenser.

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E 12. Scope of Delivery

Dispenser for soap and disinfectant

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OPHARDT HYGIENE-TECHNIK GmbH + Co. KG • Lindenau 27 • 47661 Issum • Germany